

Mixing Instructions

FORT DODGE®

Fort Dodge Animal Health

ProHeart® 6 (moxidectin)



1

Open the box and read the ProHeart 6 package insert. Remove:

- 1 vial of moxidectin microspheres (red cap) and shake to break up any aggregates
- 1 vial of vehicle for constitution (green cap)
- 25G vent needle



2

Using an 18G or 20G needle and sterile syringe, withdraw 17.0 mL of the vehicle. There is more vehicle supplied than the 17.0 mL required.



3

Insert the enclosed 25G vent needle into the microsphere vial.



4

Transfer the vehicle into the microsphere vial through the stopper using the transfer needle and syringe.



5

Remove the vent and transfer needles from the microsphere vial (properly dispose of all unused vehicle and needles). Shake vigorously until thoroughly mixed.



6

Record date and time of mixing on the microsphere vial. Let suspension stand for at least 30 minutes to allow large air bubbles to dissipate.



7

Swirl the mixture to achieve uniform suspension. Draw up the dose using an 18G or 20G needle and a 1 mL or 3 mL syringe (refer to Dosage Chart on back).



8

If administration is delayed, gently roll the dosing syringe to maintain a uniform suspension.



9

Inject subcutaneously on dorsum or side of the neck, cranial to scapula.



ProHeart[®] 6 (moxidectin)

Recommended Dosage Chart

| Dog weight (lb) | Dog weight (kg) | Dose volume (mL/dog) |
|-----------------|-----------------|----------------------|
| 11 | 5 | 0.25 |
| 22 | 10 | 0.50 |
| 33 | 15 | 0.75 |
| 44 | 20 | 1.00 |
| 55 | 25 | 1.25 |
| 66 | 30 | 1.50 |
| 77 | 35 | 1.75 |
| 88 | 40 | 2.00 |
| 99 | 45 | 2.25 |
| 110 | 50 | 2.50 |
| 121 | 55 | 2.75 |
| 132 | 60 | 3.00 |



Actual dose =
0.05 mL per kg
of body weight.